Endo Plus CL

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Thank you for choosing Endo Plus CL Precision-Tech.

All Precision-Tech brand products are guaranteed by Dentalica.
A dentist’s work is to take care of his clients, ours is to take care of every dentist.
We select the best options, we test and make them our own. Our experience guarantees their quality, both before and after the purchase.
Please do not hesitate to contact us if you have any further questions on this manual, in case of need:
Toll-free number 800.707172 - info@dentalica.com
We recommend to keep this manual in a safe place for future reference.
The manufacturer reserves the right to modify the information and the data provided in this manual when necessary and with no obligation to notify.


SYMBOLS USED IN THE USER MANUAL	
	Warning <i>If the instructions are not followed properly, operation may lead to hazards for the product or the user/patient.</i>
	Note <i>Additional information, explanation of operation and performance.</i>

SYMBOLS USED ON PACKAGING, DEVICE AND PARTS			
	Serial number		European representative
	Catalogue number		WEEE (waste of electrical and electronic equipment) directive marking
	Lot		Store in a dry place
	Manufacturer		Temperature limitation
	Date of manufacture		Relative humidity
	Consult User manual		Atmospheric pressure
	Medical Electrical Devices Directive, DS / EN 60601-1: 2006, double electrical insulation, Class II		Direct current
	Medical Electrical Devices Directive, DS / EN 60601-1: 2006 Protection of electric shocks, type B applied part		Washer-disinfector for thermal disinfection
	Medical Device Directive, MDD 93/42 / EEC, Medical Device, Class IIa		Can be autoclaved up to a maximum temperature of 134°C

1.	Scope of Endo Plus CL	pag. 42
1.1	Parts Identification	pag. 42
2.	Before Use	pag. 43
2.1	Intended Use	pag. 43
2.2	Contraindications	pag. 43
3.	Installing the Endo Plus CL	pag. 43
3.1	Installation of the contra angle	pag. 43
3.2	Install the file	pag. 44
3.3	Connecting measuring wire	pag. 45
3.4	Connecting charge base	pag. 46
4.	Use Interface	pag. 47
4.1	Panel key	pag. 47
4.2	Screen display	pag. 48
4.3	Terms and definition	pag. 50
5.	Setting	pag. 50
5.1	Selecting memory	pag. 50
5.2	Setting parameters	pag. 51
5.3	Preset programs	pag. 53
5.4	Advanced setting	pag. 54
5.5	Parameter logic	pag. 56
6.	Operation	pag. 57
6.1	Charge	pag. 57
6.2	Motor operation	pag. 58
6.3	Apex operation and not suitable condition	pag. 59
7.	Cleaning, disinfection and sterilization	pag. 63
7.1	Foreword	pag. 63
7.2	General recommendations	pag. 63
8.	Error Warning	pag. 67
9.	Troubleshooting	pag. 68
10.	Technical Data	pag. 70
10.	Components and accessories	pag. 71
12.	EMC Tables	pag. 71
13.	Statement	pag. 75

1. SCOPE OF ENDO PLUS CL

1.1 Parts Identification

1. Charge Base
2. Handpiece
3. Contra Angle
4. Insulating Sleeve
5. File Clip (2 pcs)
6. Lip Hook (2 pcs)
7. Measuring Wire
8. Adapter
9. Spray Nozzle

2. BEFORE USE

2.1 Intended Use

Endo Plus CL is a cordless endodontic treatment motorized handpiece with root canal measuring capability. It can be used to enlarge canals while monitoring the position of the file tip inside the canal. It can be used as a low-speed motorized handpiece and device for measuring canal length. This device must only be used in hospital environments, clinics or dental offices by qualified dental personnel and not used in the oxygen-rich environment.

2.2 Contraindications

The integrated apex locator of The Endo Plus CL is contraindicated in cases where patient/user carry medical implants such as pace makers or cochlear implants etc. Do not use the device for implants or other non-endodontic dental procedures. Safety and effectiveness have not been established in pregnant women and children.

Warning


Read the following warnings before use:

- The device must not be placed in humid surroundings or anywhere where it can come into contact with any type of liquids.
- Do not expose the device to direct or indirect heat sources. The device must be operated and stored in a safe environment.
- The device requires special precautions with regard to electromagnetic compatibility (EMC) and must be installed and operated in strict compliance with the EMC information. In particular, do not use the device in the vicinity of fluorescent lamps, radio transmitters, remote controls and do not use this system near the active HF Surgical Equipment in the hospital. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of The Endo Plus CL, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Do not charge, operate or store at high temperatures. Comply with the specified operating and storage conditions.
- Gloves and a rubber dam are compulsory during treatment.
- If irregularities occur in the device during treatment, switch it off. Contact the agency.
- Never open or repair the device yourself otherwise, void the warranty.


3. INSTALLING THE ENDO PLUS CL

3.1 Installation of the contra angle

Make sure 4 pins on contra angle alignment the slots of handpiece, plug them together until it “click” securely into place.


The contra angle can be 340 degrees rotated without take off, make it easy to watch the LCD in treatment by rotating the contra-angle.

Warning

- Make sure the assembly is connected properly, otherwise might cause unexpected motor reverse, even hurt the patients.
- After connecting the contra angle and handle, pull it gently to make sure the connection is good.

The improve insulation of the contra angle during combine apex, we recommend using an insulating sleeve. You can also use disposable sleeve (sold separately) instead of insulating sleeve


Note

Without the insulating sleeve, when performing the apex measurement with handpiece, wear appropriate insulated gloves, and make sure the contra angle does not touch the lips. It is advisable to use a rubber dam when performing such treatments.

3.2 Install the file

Turn the file back and forth until it is lined up with interior latch groove and slips into place, lock the file into the contra angle. Hold down the push button on the contra angle and can release the file.


Warning

- *Inspect the file head before inserting the file. Do not use the damaged file head. Make sure the motor is stopped when inserting and removing files. Be careful when inserting and removing files to avoid injury to fingers.*
- *Take care not to touch the Main switch when putting files in. this will cause the file to rotate. Pull the file gently to make sure that the file is secure in handpiece properly, otherwise it may pop out and hurt the patient.*

3.3 Connecting measuring wire

If want activity apex measurement function, uncap the USB cover on handpiece, insert measuring wire.


Insert lip hook into white slot, insert fileclip into black slot.

Note

It's not necessary to connect file clip during motor combine apex function, only during single apex function.


Note

Match colors to connect the lip hook and file clip, if connect lip hook with black slot, apex auto start will have no function.

3.4 Connecting charge base


Plug the USB of adapter into the charge base, and plug the other end into a power outlet, the Power LED on charge base will light up (green).

Note


Only the original adapter could be used.

Put the handpiece all the way into the charge base, the charge state will show on the screen.


Note

Put the handpiece into the charge base in the right direction, otherwise the handpiece will not be recharged.


If only need a base to put the device on dentist element of dental chair (without charge function), handpiece base (sold separately) is recommended to instead of charge base.

4. USE INTERFACE

4.1 Panel key

1. Main switch
2. Display Screen
3. Setting key
- 4.
- 5.

Turn Power On

Press more than 0.5 seconds to turn on the instrument.

Memory Change

Press or during standby state.

Operation mode Change

Press once during standby state, press or to change, then press or wait 5 seconds to confirm.

Parameter Adjustment

Press till target parameters, press or to adjust, then press or wait 5 seconds to confirm.

Preset Program Selection

Long press to entry preset program during standby state, press or to change, then press to confirm.


Turn Power Off

Holding down press then press .

Advanced setting


During power off state, holding down press then press to entry advanced setting, Press till target setting, press or to adjust, then press to confirm.

4.2 Screen display


Standby interface

1. Memory mode number
2. Speed of this memory program
3. Unit of speed: Revolutions Per minute
4. Residual battery power
5. Operation mode
6. Torque of this memory program
7. Unit of torque: Newton centimeter


Standard working mode

1. The set maximum torque
2. Real time torque
3. Torque display scale
4. The preset speed


Reference point interface

1. Flash bar of apical reverse position
2. Apex (Major/Anatomic apical foramen)
3. 0.5mm meter reading (Very near Minor/Physiological apical foramen)
4. 1mm-3mm (Estimate dimension) distance scale from apex (Estimate dimension)


Canal measurement standby interface

1. Memory mode number M0 is apex stand-alone memory
2. Electronic apex locator
3. Apex flash bar

Canal measurement start interface

1. Indication number. The number has no represent of actual length, only for indication
2. Canal length indicator bar
3. Apex flash bar (apex stand-alone mode) or reference point flash bar (motor combine with apex mode)

Canal measurement approach apex interface

1. The estimate dimension from apex (Major/Anatomic apical foramen)
2. Canal length indicator bar


4.3 Terms and definition

Fwd	Forward (Clockwise rotation).
Rev	Reverse (Counter clockwise rotation). Be applied to special file, inject calcium hydroxide and other solutions.
ATC	Adaptive torque control. Up to setting torque, the motor will move with reciprocating mode; when torque reduce to normal value, the motor will clockwise rotate.
EAL	Electronic apex locator. In the mode, the device will work like a stand-alone apex locator.
AP	Apex. Major apical foramen or Anatomic apical foramen.
R.L	Torque reverse less. The motor will not reverse rotation no matter how large the torque load is.
Reference point	During combined length determination, normally apical reverse must active before reaching major apical foramen, setting apical reverse position by change the flash bar.
Memory Mode	Such as M0-M10.
Operation Mode	Such as FWD, REV and ATC.


5. SETTINGS

5.1 Selecting memory

Endo Plus CL has 11 memory programs, press  or  to change during standby state, the memory number  will change according.

M1-M10 is standard working memory for canal shaping, every memory has its own speed , mode  and torque , all this these parameters can be changed (See chapter 5.2 Setting parameters).

M0 is special memory for stand-alone apex locator function (See chapter 6.3 Apex operation and not suitable condition).

5.2 Setting parameters

Warning

All the parameters must be set according to files, make sure all the parameters are expected before starting the motor, otherwise has risk of file broken.

M1 **350 rpm**
 **Fwd** **3.0 Ncm**

①

Before start the motor, check the operation mode ① is correct, otherwise press **[S]** once during standby state, press **[<]** or **[>]** to change.

Operation Mode
Fwd

Endo Plus CL has 4 operation modes: FWD, REV and ATC (See chapter 4.3 Terms and definition to get the explanations of these modes).

Note

When choice REV mode, a slow beep alarm sound appears after starting the motor, used for indicating counter clockwise rotation happening.

Repeatedly press **[S]** to check all the next level parameters of this operation mode are expected, press **[<]** or **[>]** to change if not.

Note

The parameter will differ in difference mode according to certain logic (See Chapter 5.5 Parameter logic).

Speed
300 rpm

The speed setting can be adjusted from 120 rpm to 1000 rpm.

Note

The speed of ATC operation mode is difference according to certain logic (See chapter 5.5 Parameter logic).

Torque Limit
3.0 Ncm

The torque setting can be adjusted from 0.5 N·cm to 4.0 N·cm, and R.L (torque reverse less) is also available.

Note

The torque of ATC operation mode is difference according to certain logic (See chapter 5.5 Parameter logic).

Warning

When choice R.L (torque reverse less), a slow beep alarm sound appears after starting the motor. Be careful to use this function, very professional skill is needed, otherwise has risk of file broken.

Auto Start ON

Endo Plus CL integrated apex locator, if the lip hook is connecting with patient's lip, when the endo file entering root canal, the motor will start automatically. Press or to shut off this function if not expected, press to start and stop the motor.

Warning

The motor will start automatically if handpiece (without insulating sleeve) or file touch the patient's lip or operator's fingers (without insulating glove), take care to avoid this, the file rotated by motor has risk of injure someone.

Auto Stop OFF

When the endo file out of root canal, the motor will not auto stop with default setting, Press or to select auto stop "ON" if needed.

Apical Action Reverse

Because of integrated apex locator, when the file reaches the reference point, the motor will response according to setting, it can be Reverse, SlowDown, Stop and Off.

Press or to change.

Reverse: rotation direction changing till the file upward a little bit by operator, rotation direction will change back again.

SlowDown: rotation slowdown when approach the reference point, will reverse if reach.

Stop: rotation stop when reach the reference point, upward a little bit and will rotate again.

Off: rotating as usual even if reach the reference point.


1

Reference Point

During combined length determination, normally apical reverse must active before reaching major apical foramen, Press or to set apical reverse position by change the flash bar ①, the motor will reverse while reaching the flash bar every time.


5.3 Preset programs

For convenience, we preset some common file system.

Long press **[S]** to entry preset program during standby state, the interface will show as left.

M1 **①** meanings the current memory mode, you can replace it by preset program **②** press **[<]** or **[>]** to change, then press **[•]** to confirm.


If you selecting “OneCurve” **①**, the operationmode **②**, speed **③** and torque limit **④** will change according to the file system default setting.

Note

Protaper , GATES , Pro.Glider , and Wave one is a registered trademark of Dentsply.Mtwo , Flex.Master , Reciproc and R-Pilot is a registered trademark of VDW. K3XF , TF is a registered trademark of SybronEndo. OneG , OneShape , OneFlare , 2Shape and OneCurve is a registered trademark of Micro-Mega XPendo.Shaper , XPendo.Finisher , iRace , BT-Race and BioRace is a registered trademark of FKG.


The memory mode **①** will change according, also operation mode **②**, speed **③** and torque limit **④** will charge according to the file system default setting.

Note

All of memory mode (from M1-M10) can be replaced with same method.

The parameters of “OneCurve” can also be changed make it different from default setting, and there will be 4 corners around the parameter ① ②. If want to change back to default setting, long press **[S]** to entry preset program during standby state, select “OneCurve” and press **[•]** to confirm, the default setting will be recalled, and the 4 corners around will disappear. If want to change back to M1 (or M2-M10), long press **[S]** to entry preset program during standby state, press **[<]** or **[>]** to select M1 (or M2-M10) again, then press **[•]** to confirm.

Warning

Changing the default setting is not recommended only if you are very sure it can be changed, otherwise has risk of file broken.

5.4 Advanced setting

Versions
E.1.1.008

During power off state, holding down press **[S]** then press **[•]** to entry advanced setting, the version number software will appear on the display screen. Endo Plus CL can update software very easy without tools and software. Contact your distributor to update if necessary.

Note

After updating, all of the setting parameters will be covered.

Auto Power Off
10 Min

Press **[S]** again, the “Auto Power Off” time can be change, press **[<]** or **[>]** to adjust, then press **[•]** to confirm. The “Auto Power Off” time can be set from 3-15 minutes.

Auto Return time
5 Sec

Press **[S]** again, the “Auto Return time” can be change, it means when setting parameters just like speed and torque, the system will back to standby interface if there is no operation in 5 seconds. press **[<]** or **[>]** to adjust, then press **[•]** to confirm. The “Auto Return time” can be set from 3-15 seconds.

Beeper Volume Vol.2

Press **[S]** again, the “Beeper Volume” can be change, press **[<]** or **[>]** to adjust, then press **[•]** to confirm. The “Beeper Volume” can be set from 0-3.

Habit hand Right Hand

Press **[S]** again, the “Habit hand” can be change, press **[<]** or **[>]** to adjust, then press **[•]** to confirm. The right hand and the left hand can be set.

Start up memory M1

Press **[S]** again, the “Start memory” can be change, it means every time turn power on, which memory mode will appear first. Press **[<]** or **[>]** to adjust, then press **[•]** to confirm. M1 and Last (the memory mode number when you turn power of) can be set.

Calibration OFF

Press **[S]** again, entry “Calibration” function, press **[<]** or **[>]** to select “ON”, press **[•]** to start calibration.

Warning

Before calibrating, making sure the original contra angle is installed, and do not install the file. The torque will not correct if calibration without original contra angle or any load on contra angle chuck, and has risk of file broken.

Calibration 1000 rpm

The motor speed will increase from 120 to 1000 rpm. When the speed up to 1000 rpm, the calibration successful and automatic power off.

Restore settings OFF

Press **[S]** again, entry “Restore setting” function, press **[<]** or **[>]** to select “ON”, press **[•]** to start recovering, all the parameters be set by operator will be recovered by default factory setting (See chapter 5.5 Parameter logic).

Note

After restore setting, all the parameters will be covered, record what you need before this operation.

5.5 Parameter logic

The default memory settings are listed below, the setting can be changed as needed.

Function	M1	M2	M3	M4	M5	M6	M7	M8	M9	M10
Operation Mode	FWD	FWD	FWD	FWD	REV	REV	REV	REV	FWD	FWD
Speed (rpm)	300	400	500	300	300	300	350	500	800	1000
Torque Limit (N·cm)	3.0	2.0	2.0	3.0	1.5	1.0	2.5	2.0	1.5	1.0
Auto Start	ON	ON	ON	ON	ON	ON	ON	ON	ON	ON
Auto Stop	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF
Apical Action	REV	REV	REV	REV	REV	REV	REV	REV	REV	REV
Reference point	02	02	02	02	02	02	02	02	02	02
FWD Angle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
REV Angle	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

The default advanced settings are listed below, the setting can be changed as needed.

Auto Power off	10 Min
Auto Return time	5 Sec
Beeper Volume	2
Habit hand	Right hand

Startup memory	M1
Language	English
Calibration	OFF
Restore settingsw	OFF

The speeds (rpm) in different operation mode are not the same, details are listed below.

FWD	REV
120 150 200 250 280 300 350 400 450 500 550 600 650 700 750 800 850 900 950 1000	

The torques (N·cm) in different operation mode are not the same, and even in the same operation mode, when the speed changing, the possible torque is difference, details are listed below.

Speed (rpm)	FWD	REV									
120-700	0.5 0.8 1.0 1.5 1.8 2.0 2.2 2.5 3.0 3.2 3.5 4.0 R.L										
700-1000	0.5 0.8 1.0 1.5 1.8 2.0										

6. OPERATION

6.1 Charge

Displays the present remaining amount of the battery. Less than 15% remains, please charge.


Note
If the power is less than 15%, must be recharged within 30 days, otherwise the battery will be damaged.

Low Power
Please Charge

If continue to use, the torque and speed will low than the setting value, and low power warning will appear on screen, and the device will stop work.

Note
The remaining amount of battery mark indicates a voltage. When a load is applied to the motor handpiece, the remaining amount of battery mark appears to become lower.


Alternative charging method

Charge without charge base also available, using adapter connect to handpiece directly, the charge state will show on the screen. Charge with charge base is recommended (See chapter 3.4 Connecting charge base).

Note
Only the original adapter could be used.

Charging indication appears on the screen, and flashes slowly ①, when battery is fully charged or in a state near full charge, the flash will stop and show like picture ②. Fully charged will take about 4 hours, depending on residual battery power and battery state. It can be recharged 300-500 times, depending on the operating conditions of the device.

Note


When changing, other function will forcibly stop, take from charge base, press main switch, the last function will recall.

Warning

Do not change the battery, only trained technician or distributor can change the battery, the electronic parts will be damaged if use a wrong battery or install with a wrong way.

6.2 Motor operation

When using as a stand-alone motor, the torque bar will show on the screen (more information about torque bar, please see chapter 4.2 Screen display).

Warning

- Use The Endo Plus CL outside the oral cavity to make sure that the device is functioning properly.
- Change file on time to avoid file separation within the canal. File may separate because of cyclic / torsional fatigue.
- Heavy force / hand pressure on endo motor while using may even cause file separation.
- Do not press the button to release the files while the motor is running, otherwise the file may pop out and even hurt the patient.
- Electromagnetic noise in surroundings environment may interfere with the device operation, do not rely on device's automatic control completely, always pay attention to the feedback from display.



Note

If there is any abnormal functioning, stop using the device and report to company. The file separates more easily at high speeds, please follow the manufacturer's recommendations of the speed and check the settings of the Endo motor before use. Do not use the files are except nickel-titanium or stainless steel. Gloves and a rubber dam are compulsory during treatment. Do not forget to remove the file from the Contra-angle after its use.

When using motor combine apex function, the measure wire must be connecting with motor by USB socket, and white slot connects with lip by lip hook, keep the black slot idle. The reference point bar will show on the screen (*more information about reference point bar, please see chapter 4.2 Screen display*).

Note


We strongly recommend check the function every time before use. Touch the lip hook with the file in the contra angle and check that all the bars on the meter on the screen light up, and the motor should be reversed continuously.


Note

The will not be able to perform a precise measurement for every time, especially in cases of abnormal or unusual morphology of the root canal. The user need coordinate with x-ray to check the results of the measurement. If the meter does not move when you enter the file, it is possible that the unit is not working normally, therefore, stop using.


6.3 Apex operation and not suitable condition

When using as a stand-alone apex locator. We suggest put the handpiece on the charge base to get better visual angle. The measure wire must be connecting with motor by USB socket, white slot connects with lip by lip hook, and black slot connect with file clip. the canal length indicator bar will show on the screen (*more information about canal length indicator bar, please see chapter 4.2 Screen display*).


The reference point can be adjusting when use stand-alone apex function. Press **[S]** to active reference point interface during M0 standby state, Press **[<]** or **[>]** to change reference point by change the flash bar **①**, a continuous beeping with appear when reach the reference point.


Note

The file clip must hold the file correctly. Push the button of file clip in the direction shown by the arrow, clip the holder onto the metal upper part of the file and then release the button. The file clip must be nearly vertical with file handle, otherwise will damage the tip of the file holder.

Note

We recommend using Apex Tester to check the accuracy of apex locator every time before use. Insert the Apex Tester into handpiece USB socket during M0 mode, check the highlight number must between 01-03 (0.3mm-0.8mm on the top).


Note

We recommend check the connection of apex locator every time before use. Touch the lip hook with metal part of file clip, check all the bars on the meter on the screen light up, and "over" flash on the top.


Unsuitable situation of root canals for Electric Measurement

Cannot obtain precise measurements if the root canal conditions as below

Root canal with a large apical foramen

The root canal cannot be accurately measured because of the lesion or incomplete development of the apical foramen. The results may show that the length measured is shorter than the actual one.


Root canal blood overflow from the opening

If blood spills from the root opening and contacts the gums, it will cause leakage of electricity, which cannot be accurately measured. Wait for the bleeding to stop completely. Clean the root canal and the opening, completely empty the root canal blood, and then measure it.

The root canal uses a chemical solution to flow out from the opening

If a chemical solution flows out of the root canal, it is impossible to get an accurate measurement. It is important to remove the overflow from the opening.


Build-up (e.g. cement)

Broken crown

If the crown is broken, a segment of the gingival tissue enters the lumen, and the contact between the gingival tissue and the root file causes electrical leakage, which cannot be accurately measured. In this case, the appropriate material should be used to isolate the gingival tissue.


Crack

The crack tooth Leakage through branch of the root canal

Broken teeth can cause electrical leakage and cannot be accurately measured. Branch tubes can also cause leakage.


Branch

Retreatment canal which was filled with gutta-percha


The gutta-percha must be completely removed to eliminate its insulation, then pass a small file all the way through the apical foramen and then put a little saline in the canal, but do not let it overflow the canal opening.

Crown


Crown or metal prosthesis that touches gingival tissue

Accurate measurement cannot be obtained if the file touches a metal prosthesis that is touching gingival tissue. In this case, widen the opening at the top of the crown so that the file will not touch the metal prosthesis before taking a measurement.


Cutting debris on tooth Pulp inside canal

Remove all cutting debris on the tooth. Remove all the pulp inside the canal. Otherwise an accurate measurement cannot be obtained.


Caries touching the gums

In this case, electrical leakage through the caries infected area to the gums are impossible to obtain an accurate measurement.

Blocked canal

The meter will not run if the canal is blocked. Opening the canal all the way to the apical construction to measure it.


Extremely dry canal

If the canal is extremely dry, the meter may not work until it is quite close to the apex. In this case, try to moisten the canal with oxydol or saline.

Difference measuring result between Apex locator reading and Radiography

Sometimes the reading of the apex locator reading does not correspond to the Xray image. this does not mean inaccurate of apex locator or X-ray, depending on the angle of the X-ray beam, the root tip may not be displayed correctly. The position of the root tip seems to differ from its true position.

The X-ray photo shows that the actual apex of the root canal is not the same as the anatomic end. In fact, the apical foramen is located at the coronal end. in this case, X-ray may indicate that the file needle has not reached the apical foramen, even if it has actually reached the apical foramen.

7. CLEANING, DISINFECTION AND STERILIZATION




7.1 Foreword

For hygiene and sanitary safety purpose, the components (contra angle, file clip, lip hook and insulating sleeve) must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent uses.

Comply with your national guidelines, standards and requirements for cleaning, disinfection and sterilization.

7.2 General recommendations

- 7.2.1 The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments, where applicable after sterility.
- 7.2.2 For your own safety, please wear personal protective equipment (gloves, safety glasses, etc.).
- 7.2.3 Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, and FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer.
- 7.2.4 The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- 7.2.5 Thoroughly clean and wash the components before autoclaving.
- 7.2.6 Do not lubricate the motor handpiece.
- 7.2.7 Do not clean the contra angle with an ultrasonic cleaning device.
- 7.2.8 Do not use bleach or chloride disinfectant materials.


AUTOCLAVABLE COMPONENTS			
CONTRA ANGLE	LIP HOOK	FILE CLIP	INSULATING SLEEVE
			



Warning

- Only the components above can be autoclaved.
- Before first use and after each use, sterilize the above components.

STEP NO.	INSTRUCTIONS	
1	Initial treatment at point of use	Immediately after using, wipe gross contaminations from the components (contra angle, file clip, lip hook and insulating sleeve), and put them in container for transportation. Prepare the components directly after treatment.  Warning Do not submerge the components or wipe them with any of the following functional water (acidic electrolyzed water, strong alkaline solution, or ozone water), medical agents (glutaral, etc.), or any other special types of water or commercial cleaning liquids. Such liquids may result in metal corrosion and adhesion of the residual medical agents to the components.
2	Preparation before cleaning	Remove and disconnect the components (Contra Angle, Lip Hook, File clip, Insulating Sleeve) before cleaning (refer to “Chapter 3-Installing the Endo Plus CL” of this manual for disassembly instructions).  Warning - Do not fail to take out the file before cleaning the contra angle. - Observe suitable personal protective measures.
The following Step 3 to Step 5 are operated in a washer-disinfector:  Warning - Use only approved washer-disinfectors according to EN ISO 15883, maintain and calibrate it regularly. - Follow instructions and observe concentrations given by the manufacturer (see general recommendations). - Sufficient rinsing step should be available in purified water (max 10 germs/ml and max 0.25 endotoxin units/ml). - Avoid any contact between the contra-angle and any instrument, kit, support or container. - Make sure the components are dry before moving to the #6 step.		

STEP NO.	INSTRUCTIONS	
3	Cleaning: Automated	• Carefully put the components (Contra Angle, Lip Hook, File Clip, Insulating Sleeve) into the washer-disinfector and set the parameters as follows: • Pre-cleaning: water temperature <30°C, 2 min; • Cleaning: water temperature 45°C, 5 min; use an enzyme detergent solution (mild and aldehyde free solution) which is suitable to be used with washer-disinfector, and use in accordance with the IFU of the detergent solution manufacturer; • Rinsing: water temperature 45°C, 1 min (Rinsing twice).
4	Disinfection: Thermal	• Thermal disinfection at least 5 min at 90°C/194°F, make sure A0 value ≥ 3000.
5	Drying	Heat, 20 min, 90°C/194°F.
6	Maintenance and Inspection	• Inspect components and sort out those with defects. Dirty components must be cleaned and disinfected again. • Only the contra angle needs to be lubricated.  Black oil  Warning • Before autoclaving, the contra angle must be lubricated. • Attaching the spray nozzle to oil can and contra angle, press the oil can button more than 3 seconds, till all the black oil flow out from the head of the contra angle.
7	Packaging	Pack each component in a separate steam-sterilization pouch.  Warning • Check the validity period of pouch given by the manufacturer to determine the shelf life. • Use pouches which resist to a temperature up to 141°C (286°F) and in accordance with EN 868.

8	Sterilization	• Steam sterilization at 134°C, at least 5 minutes. • Minimum drying time after sterilization: 10 minutes.  Warning • Use only approved autoclave devices according to EN 13060 or EN 285. • Use a validated sterilization procedure according to ISO 17665. • Respect the maintenance procedure of the autoclave device given by the manufacturer. • Use only this recommended sterilization procedure. • Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physicochemical integrators, digital records of cycles parameters). • The sterilization procedure must comply with ISO 17665. • Wait for cooling before touching.
9	Storage	Keep the components in sterilization packaging in a dry and clean environment.  Warning • Sterility cannot be guaranteed if packaging is open, damaged or wet. • Check the packaging and the contra angle before using it (packaging integrity, no humidity and validity period).
 Note The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.		

DISINFECTION COMPONENTS		
MOTOR HANDPIECE	CHARGE BASE	ADAPTER
		
MEASURING WIRE		
		

Wipe all the surfaces with a cloth lightly moistened with Ethanol for Disinfection (Ethanol 70 to 80vol%) at least 2min, repeat for 5 times.

Note

Do not use anything except Ethanol for Disinfection (Ethanol 70 to 80 vol%). Do not use too much ethanol as it's going into machine and damage the components inside.

8. ERROR WARNING

Overload Restart Motor

Overload

When setting the torque limit as R.L or during reverse processing, the Overload warning may appear on the screen, it means a large load happened greater than the motor force. Press the Main switch to restart motor.

Overheat See user manual

Overheat

The temperature of motor is higher than expectation, turn the power off and waiting more than 5 minutes to let it cold down.

HWFault See user manual

HWFault

Hardware of the handpiece broken, contact your distributor.

Motor Fault
See user manual

MotorFault

Motor of the handpiece broken, contact your distributor.

LowPower
Please Charge

LowPower

The power is very low, charge it immediately

9. TROUBLESHOOTING

When trouble is found, check the following points before contacting your distributor. If none of these are applicable or the trouble is not remedied even after action has been taken, the product may have failed. Contact your distributor.

PROBLEM	CAUSE	SOLUTION	REF. CHAP
The power is not turned on.	The battery is flat.	Charge the battery.	6.1
	Press  too short time.	Press  more than 0.5 seconds.	4.1
The power led on charge base does not light.	Using a wrong adapter.	Use the original adapter.	3.4
	The adapter is not connected.	Check the connection.	3.4
	The plug of the adapter is not inserted into the outlet.	Check the connection.	/
	There is no electricity in the outlet.	Check the connection.	/
No charge indicator flash on handpiece screen	Put the handpiece into the charge base in the wrong direction.	Check the direction.	3.4
	Charge pin of charge base unable to rebound.	Remove debris which between move part and base of the charge pin.	/
	Contactors are dirty.	Cleaning the surface of contactors.	/
	The charge base broken.	Using adapter connect to handpiece directly, and Contact your distributor.	/
Handpiece screen does not appear	The handpiece broken.	Check if there is a sound of beep or motor, and Contact your distributor	/





PROBLEM	CAUSE	SOLUTION	REF. CHAP
The motor doesn't rotate	M0 mode is stand-alone apex locator function.	Changing to M1-M10.	5.1
	The contra-angle is clogged.	Clean or replace the contra-angle.	/
	Motor is protected by system or broken.	Check the error warning.	8
Motor does not run when the file is inserted in the canal.	The measure wire connecting not properly.	Check the connection.	3.3
	The lip hook not properly hooked in the corner of the patient's mouth.	Check the connection.	6.2 6.3
	The AUTO START function is OFF.	Turn the auto start function ON if necessary.	5.2
The motor Can't stop	The AUTO STOP function is OFF.	Press  to stop it, setting AUTO STOP function ON if necessary.	5.2
	There is a short circuit inside the motor handpiece or the motor handpiece cord.	Press  button to stop the motor and contact your distributor.	/
Motor spontaneously starts running in reverse.	Up to setting torque limit.	Check the torque limit is enough or not.	5.2
	Apical action setting to Reverse	Change setting if it's not expected.	5.2
	Setting to REV mode.	Change setting if it's not expected.	5.2
Motor does not reverse.	R.L mode is set.	Change setting if it's not expected.	5.2
	Torque reverse setting might be too high.	Change setting if it's not expected.	5.2
	Apical action setting Stop or OFF.	Change setting if it's not expected.	5.2
Motor speed changes spontaneously.	Apical action setting Slow Down.	Change setting if it's not expected.	5.2
Motor alternates between forward and reverse rotation.	Operation mode setting to ATC.	Change setting if it's not expected.	5.2
No sound.	Beep volume set to 0.	Set beep volume to 1, 2 or 3.	5.4
Beep sound an alarm even though the instrument is not being used.	The motor is set to REV or R.L mode.	If it is the expected mode, ignore the alarm.	5.2

PROBLEM	CAUSE	SOLUTION	REF. CHAP
Canal measurement is unstable.	Complex root canal environment.	Check situation of root canals.	6.3
Cannot make a measurement.	Measure wire, lip hook or file clip connecting not properly.	Check the connection.	6.2 6.3
	Lack electrical conductivity between the shank and the file.	Use a file that has conductivity.	/
	Unsuitable situation of root canals.	Check the root canal environment.	6.3


10. TECHNICAL DATA

MANUFACTURER	Changzhou Sifary Medical Technology Co., Ltd
MODEL	Endo Plus CL
DIMENSIONS	21.5cm x 17.5cm x 9cm±1cm (Outer box)
WEIGHT	1 kg ± 10%
CONTRA-ANGLE	Contra-angle compatible with rotary and reciprocating instruments, equipped with a 2.35 mm shaft conforming to ISO 1797-1:2011, Type 1, Files length 11-31mm.
POWER SUPPLY	Lithium ion battery: 3.7V, 1500mAh, ±10%
CHARGER POWER SUPPLY	AC 100-240 V, ± 10%
CHARGE POWER OUTPUT	5V $\pm$ 1A
FREQUENCY	50/60 Hz, ± 10%
CHARGER NOMINAL POWER INPUT	5.5VA
TORQUE RANGE	0.5 Ncm – 4 Ncm
SPEED RANGE	120-1000 rpm
TYPE OF PROTECTION AGAINST ELECTRICAL SHOCK	Class II and internally powered equipment
APPLIED PART	B
OPERATION CONDITIONS	Use: in enclosed spaces Ambient temperature: 5°C - 40 °C Relative humidity: <80% Operating altitude < 3000m above sea level Atmospheric pressure: 70kPa - 106kPa
TRANSPORT AND STORAGE CONDITIONS	Ambient temperature: -20 °C - +55 °C Relative humidity: 20% - 80 % Atmospheric pressure: 70kPa - 106kPa

11. COMPONENTS AND ACCESSORIES

MOTOR HANDPIECE (1PCS)	CHARGE BASE (1PCS)	CONTRA ANGLE (1PCS)
		
	PR00ENDOPCB	PR00ENDOPCT
MEASURING WIRE (1PCS)	LIP HOOK (2PCS)	FILE CLIP (2PCS)
		
PR00ENDOPCCV	PR00ENDOPCG	PR00ENDOPCC
INSULATING SLEEVE (1PCS)	SPRAY NOZZLE (1PCS)	ADAPTERS
		
PR00ENDOPCCCI		

OPTIONS (SOLD SEPARATELY)

HANDPIECE BASE	APEX TESTER (1PCS)
	
PR00EPCGPPM	PR00ENDOPCTR

12. EMC TABLES

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
RF emissions CISPR 11	Group 1	The Endo Plus CL uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
IEC 61000-3-3	Complies	
		The Endo Plus CL is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Endo Plus CL is intended for use in the electromagnetic environment specified below.
The customer or the user of The Endo Plus CL should assure that it is used in such an environment.

TEST DI IMMUNITÀ	IEC 60601 TEST	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	+/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transients/bursts IEC 61000-4-4	±2kV 100kHz repetition frequency	±2kV 100kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Line to line: ±0.5kV, ±1kV Line to earth: ±0.5kV, ±1kV, ±2kV	Line to line: ±0.5kV, ±1kV Line to earth: ±0.5kV, ±1kV, ±2kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips IEC 61000-4-11	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0°	0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles sine phase at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of devices require continued operation during power mains interruptions, it is recommended that devices be powered from an uninterruptible power supply or a battery
Voltage interruptions IEC 61000-4-11	0% UT; 250/300 cycle	0% UT; 250/300 cycle	
Rated Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz or 60Hz	30 A/m 50 Hz o 60 Hz	Power frequency magnetic field should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz or 30 cycles at 60Hz

The Endo Plus CL is intended for use in the electromagnetic environment specified below.
 The customer or the user of The Endo Plus CL should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST	LIVELLO DI CONFORMITÀ	AMBIENTE ELETTROMAGNETICO GUIDA
Conducted dis-turbances induced by RF fields IEC 61000-4-6	3 V 0.15 MHz – 80 MHz, 6 V in ISM bands between 0.15 MHz and 80 MHz, 80 % AM at 1 kHz	3 V	Portable and mobile RF communications equipment should be usedno closer to any part of The Endo Plus CL, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF EM fields IEC 61000-4-3	3 V/m, 80 MHz – 2,7 GHz, 80 % AM at 1 kHz	3 V/m	
Proximity fields from RF wireless communication equipment IEC 61000-4-3	See the RF wireless communication equipment table in “Recommended minimum separation distances”	Complies	Recommended minimum separation distances. See the RF wireless communication equipment table in “Recommended minimum separation distances”

RECOMMENDED MINIMUM SEPARATION DISTANCES

Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The Endo Plus CL has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipments and The Endo Plus CL as recommended below.

TEST FREQUENCY (MHz)	BAND (MHz)	SERVICE	MODULATION	MAXIMUM POWER (W)	DISTANCE (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430-470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 - 745- 780	704-787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
810 - 810 - 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
1720 - 1845 - 1970	1700- 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
2450	2400- 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240, 5500, 5785	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9

Warning

- Use of accessories and cables other than those specified or provided by the manufacturer of Endo Plus CL could result in increased electromagnetic emissions or decreased electromagnetic immunity of Endo Plus CL and result in improper operation. Cable information:

Cable Name	Cable Length (m)	Shielded or not	Remark
Adapter Cable	1.2	No	/
Measuring Wire	1.5	No	/

- Use of Endo Plus CL adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, E-connect S and the other equipment should be observed to verify that they are operating normally.

13. STATEMENT

Disposal


The package should be recycled. Metal parts of the device are disposed as scrap metal. Synthetic materials, electrical components, and printed circuit boards are disposed as electrical scrap. The lithium batteries are disposed as special refuse. Please deal with them according to the local environmental protection laws and regulation.

Rights

All rights of modifying the product are reserved to the manufacturer without further notice. The pictures are only for reference. The final interpretation rights belong to CHANGZHOU SIFARY MEDICAL TECHNOLOGY CO., LTD. The industrial design, inner structure, etc, have claimed for several patents by SIFARY, any copy or fake product must take legal responsibilities.

Warranty

Endo Plus CL is covered by 1 year warranty.

Selezionato, distribuito e garantito da
DENTALICA S.p.A. Via Rimini 22 | 20142 Milano - Italia
Tel. +39 02 89.598.1 | Fax +39 02 89.50.42.49
www.dentalica.com | dentalica@dentalica.com

Changzhou Sifary Medcial Technology Co., Ltd
Add: NO.99, Qingyang Road, Xuejia County, Xinbei
District, Changzhou City, 213000 Jiangsu, China
Tel: +86-0519-85962691
Fax: +86-0519-85962691
Email: ivy@sifary.com
Web: www.sifary.com

Llins Service & Consulting GmbH
Tel: +49 175 4870819
Add: Obere Seegasse 34/2, 69124,
Heidelberg, Germany
Email: Llins.Service@gmail.com