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Composite for Bulk Fill technique

INSTRUCTIONS FOR USE

General
Bulk-Fill Composite is a light curing, flowable, highly radiopaque (210 % Al), composite for posterior restorations, for direct filling and the bulk fill technique. The product complies with the guideline DIN EN ISO 4049.

Composition
Monomer matrix: diurethane dimethacrylate, butanediol dimethacrylate Total Filler: 77 % (by weight) glass filler, pyrogenic silicic acid.

Indications
- Fillings with layer thicknesses up to 4 mm in Black's class I and II cavities.
- Cavity lining - as a first (bottom) layer in Black's class I and III cavities.

Contraindications/Interactions
In case of hypersensitivities to a component, it is recommended to not use this product or to use it under strict medical supervision by the doctor/dentist. Before use, dentist has to consider all the interactions and cross-reactions of this product with other materials already present in the patient's mouth. Unpolymerized composite may cause skin sensitization in allergic or susceptible individuals. The user has to take adequate precautions (such as gloves) and be informed about the case history of the patient. In case of irritation or allergy due to one of the constituents listened, do not use this product.

Side effects
Unwanted side effects are extremely rare if used in a proper manner. To prevent possible reactions of the pulp into cavities where the dentine is exposed, the pulp must be adequately protected (such as calcium hydroxide preparation).

Application
Preparatory measures
Before commencing the treatment, clean the tooth with non-fluoride polishing paste.

1. Cavity preparation
Prepare the cavity minimally invasively as generally required for adhesive techniques. All enamel margins in the anterior region must

Per mayor información dirigirse a:
DENTALICA S.p.A. – Via Rimini 22 – I-20142 Milano – Italia
Tel.#39 - 02895981 – Fax #39 - 0289504249
www.dentalica.com – E-mail: dentalica@dentalica.com

Per ulteriori informazioni rivolgetevi a:
DENTALICA S.p.A. – Via Rimini 22 – I-20142 Milano – Italia
Tel.#39 - 02895981 – Fax #39 - 0289504249
www.dentalica.com – E-mail: dentalica@dentalica.com

Rev 0_07/13

be beveled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.

2. Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

3. Etching with etching gel

Ensure a very precise application of etching gel (for ex. Axia Etch Precision), beginning with the enamel margins. Allow it to react for at least 15 seconds before spraying thoroughly with water and drying the etched surfaces with oil-free compressed air. Avoid dehydrating the dentine. After drying, the etched surfaces of the enamel appear chalky/white and must not be contaminated prior to applying the bonding agent. Should these surfaces be contaminated with saliva, rinse and dry them again, re-etch if necessary.

4. Applying the bonding agent

Refer to the adhesive's instructions for use (for ex. One-Q-Bond C.G.T. Precision).

5. Application of the composite

Place layers of max. 4 mm of A-Bulk composite directly into the cavity using the curved application tips supplied with the material. For hygiene reasons these application tips are must only be used once and then discarded. When applying the material with the tip, make certain that no air bubbles become entrapped in the material. Wet the prepared tooth surfaces thoroughly. Polymerize each layer with a conventional polymerization unit for 40 seconds holding the light guide as close as possible to the surface. During polymerization, a dispersion layer forms on the material surface. This layer must not be touched or removed. It forms the bond with the layer of composite applied subsequently (e.g. with a universal or posterior composite) or for creating aesthetic cusps. Alternatively, A-Bulk can be applied in a thin layer to form a liner on exposed dentin. The material can be polymerized in layers of up to 4 mm thickness.

6. Finishing

A-Bulk is finished and polished with finishing diamond tools, flexible discs, silicone polishers and polishing brushes immediately after polymerization.

Notes

When creating large and time-consuming restorations either cover the material with a light proof foil or temporarily move the operating light away from the site in order to avoid premature polymerization of the material.

For hygienic reasons the curved application tips supplied with the material must be used only once

Use a polymerization unit with an emission spectrum of 350-500 nm to cure this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the manufacturer.

Precautions

Close the syringes tightly immediately after use. After use retract the plunger of the syringes lightly to prevent the apertures from becoming blocked. Let the material reach room temperature before use. This product was developed specifically for the described range of application. It must only be used in accordance with the above instructions. The manufacture will not be held liable for damages caused by incorrect handling or incorrect processing of the material.

Storage

Store the product < 25°C (77°F). Avoid direct sunlight. Do not use after expiration date.

Packaging

2 syringes of 2 gr each, 20 application tips.

For dental use only. Keep out of reach of children.

PRECISION products are carefully developed and extensively tested. We therefore guarantee their impeccable quality when used in accordance with their instructions for use. Application, use and processing of the product are beyond our control and are subject to the responsibility of the user.

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