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ONE-Q-BOND C.G.T.

Universal, no-rinse, self-priming etchant adhesive system

EN INDICATIONS

1. Etch dentin for 5-10 seconds. For enamel, etch for 20-30 seconds.
2. Rinse with a spray of oil-free water for 5-10 seconds.
3. Blot dry excess moisture from the surface. Verify that the surface has no visible moisture.
4. Dispense One-Q-Bond C.G.T. directly onto a microbrush
5. Apply ample amount of One-Q-Bond C.G.T. directly onto a microbrush.
6. Remove solvents with a gentle air flow for at least 15 seconds or until no visible movement of the resin is noted.

If the surface appears incompletely partially coated, apply additional material and dry again.

7. Light-cure One-Q-Bond C.G.T. for 20 seconds. If the restoration retention area lies primarily within dentin, apply a second layer of One-Q-Bond C.G.T., remove solvents immediately and light-cure for 20 seconds.
8. Restore with light-cure or self-cure composite resin of your choice or cement restoration with Composite resin cement. For light-inaccessible procedures, One-Q-Bond C.G.T. can be self-cured by utilizing the One-Q-Bond Activator. Mix equal number of drops of One-Q-Bond C.G.T. and Activator in a clean mixing well for 2 seconds and proceed with step 5.

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QUASER S.r.l.
Via Pian due Torri 19, 00146 Rome, Italy

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