ONE-Q-BOND C.G.T.

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ONE-Q-BOND CGT

Universal, no-rinse, self-priming etchant adhesive system

INDICATIONS

- Direct composite restorations.
- Adhesive cementation of crowns and bridges, inlays and onlays.

Properties

- One-Q-Bond C.G.T. is a 6th generation self-etch bonding system, which delivers high bonding values to both dentin and enamel.
- One-Q-Bond C.G.T. consists of an acidic water based self-etch Primer and a light-cure Adhesive.
- It is used to bond composite resin filling materials, to tooth structure, any self-cure or light-cure composite, compomer restorative materials and composite resin cements available on the market.
- One-Q-Bond C.G.T. is utilizing a moderately acidic Primer. Unlike total-etch, total-bond techniques, the self-etch Primer does not open the dentinal tubules completely. The smear layer is solubilized, and due to the high hydrophilic properties of the Primer, it is then able to penetrate to the tubules and peritubular dentin, forming resin tags. In the case of enamel, the Primer produces a significant pattern with enhanced surface area, leading to improved enamel bonding.
- One-Q-Bond C.G.T. has been on the market since 2008 with no changes in formula. No complaints of any health-related or chemical performance issues were received.

Caution!

- One-Q-Bond C.G.T. contains polymerizable monomers which may cause allergic contact dermatitis in susceptible individuals. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if you know allergy to methacrylate resin exists, discontinue use.
- Polymerization may be triggered by ambient lighting or the dental operator light. Tightly recap the bottle after use.
- One-Q-Bond C.G.T. contains materials which: sensitive to light/ highly flammable / irritating to eyes / harmful if swallowed / irritating to respiratory system and skin / may cause sensitization by skin contact / repeated exposure may cause skin dryness or cracking / vapors may cause drowsiness and dizziness. For complete symbols guide see below.

Composition

Essential ingredients - ETCHING	% Wt.	Essential ingredients - BOND	% Wt.
Multifunctional monomers	<20%	Multifunctional methacrylates	<30%
Methacrylated phosphoric acid esters	<5%	Multifunctional monomers	<10%
Adhesive accelerator	<10%	Urethane dimethacrylate	<10%
Water	<40%	Photoinitiators	<1%
Ethanol	<30%	Acetone	<50%
Acetone	<10%	Ethanol	<10%
		Adhesive accelerator	<5%

Filler content

The product does not contain filler.

INSTRUCTIONS FOR USE

Application of One-Q-Bond C.G.T. Self-Etching:

1. Apply ample amount of One-Q-Bond C.G.T. SE on the prepared, air-dried or blotted tooth surface using a sponge or disposable brush tip. Leave undisturbed for 15 seconds.
2. Remove solvents with a blast of air for at least 10 seconds.
3. Tightly re-cap bottle after use.

Application of One-Q-Bond C.G.T.:

1. Apply ample amount of One-Q-Bond C.G.T. using a sponge or disposable brush tip.
2. Remove solvents with a gentle air flow for at least 10 seconds or until no visible movement of the resin is noted. If the surface appears partially coated, apply additional material.
3. Tightly re-cap bottle after use.
4. Light-cure for 20 seconds with a light cure unit.
5. When the restoration retention area lies primarily within dentin, apply a second layer of One-Q-Bond C.G.T. Remove solvents with a gentle air flow for at least 10 seconds. Do not cure.
6. When light curing is not an option (bonding in root canal), add equal drops of One-Q-Bond Activator with One-Q-Bond C.G.T. in a clean dappen dish mix for 2 seconds and proceed with steps 1, 2.
7. Restore with light-cured or self-cure composite resin of your choice or cement restoration.

Storage and disposal

- Do not store above 24°C (75°F).
- Keep product in original packaging.
- Keep away from direct sunlight.
- For optimum freshness, keep refrigerated.
- Do not freeze.
- To obtain SDS or IFU visit www.dentalica.com.
- Do not use after expiration date. See expiry date on the label.
- Shelf life is 2 years after the date of manufacturing.
- Empty product can be safely disposed after use in an intended regular waste container. The product does not require special disposal instructions or precautions and is harmless to the environment.

Packaging

One-Q-Bond C.G.T.: 2 Bottles of 5 ml

One-Q-Bond SE: 2 Bottles of 5 ml

For professional dental use only. Made in Israel.

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MODE D'EMPLOI

Application de One-Q-Bond C.G.T. Self-Etching (SE):

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B.J.M. Laboratories Ltd.
Dentalica S.p.A., Via Rimini 22, I-20142 Milano, Italia