ONE-Q-BOND C.G.T. Self Etching Agent

Primer idrofilo auto-mordenzante

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ONE-Q-BOND C.G.T. Self Etching Agent

Universal, no-rinse, self-priming etchant adhesive system

EN INDICATIONS

- Direct composite restorations.
- Adhesive cementation of crowns and bridges, inlays and onlays.

1. Apply ample amount of One-Q-Bond SE on the prepared, airdried or blotted tooth surface using a sponge or disposable brush tip. Leave undisturbed for 15 seconds.
2. Remove solvents with a blast of air for at least 10 seconds. Tightly re-cap bottle after use.
3. Apply ample amount of One-Q-Bond C.G.T. Bond using a sponge or disposable brush tip
4. Remove solvents with a gentle air flow for at least 10 seconds or until no visible movement of the resin is noted. If the surface appears partially coated, apply additional material. Tightly recap bottle after use.
5. Cure for 20 seconds with a visible light cure unit.

6. When the restoration retention area lies primarily within dentin, apply a second layer of One-Q-Bond C.G.T. Bond. Remove solvents with a gentle air flow for at least 10 seconds. Do not cure.
7. When light curing is not an option (bonding in root canal), add equal drops of One-Q-Bond Activator with One-Q-Bond in a clean dappen dish mix for 2 seconds and proceed with steps 1, 2.
8. Restore with light-cured or self-cure composite resin of your choice or cement restoration.

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