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DETAIL MONO NEXT

Detail curing silicone for monophase impression with high stability and fluidity

INSTRUCTIONS FOR USE

General

Detail Mono Next is a precision impression material with medium consistency for mixing devices, based on Vinyl polysiloxanes (VPS) for addition.

Indications

- Impressions on natural teeth for bridges and crowns, inlay and onlay
- Functional and relin impressions
- Implant impressions with pick-up techniques or transfer relocation

Remarks

None known to date. It cannot be excluded that especially predisposed persons might, in singular cases, develop hypersensitivity against single components of the material.

Do not use latex gloves. Avoid contact of the material with organic solvents or liquids containing organic solvents, this could modify its characteristics and cause

a lack of precision. Moreover, possible residues of retraction aids can increase the setting time. Do not use the material with condensation curing silicone.

Impression trays and adhesive

Detail Mono Next can be used with non-perforated, metal or resin trays or individual trays. Be sure using a tray with adequate space available between the tray wall and dental arch. It is possible to use a tray adhesive for VPS materials in a thin and uniform film, in order to improve the adhesion between Detail Mono Next and nonperforated tray. Let the adhesive dry at least 3 minutes on the tray before inserting the impression material.

Disinfection and pouring the impression

The impression can be disinfected with commercially available immersion solutions. We recommend using a 2% glutaraldehyde solution. Ensure that the disinfectant is indicated for use with additioncured impression silicones.

Immerse the impression for approx. 10 to 15 minutes. Then rinse under running water for 15 seconds. Keep the impression at normal room temperature (max. 25°C) and pour the impression no sooner than 30 minutes after removing it from the patient's mouth.

The impression material is compatible with all commercial casts. We recommend using a special purpose hard plaster (ISO 6873, Type 3; e.g. Moldano, Bayer) or an ultra hard dental plaster (ISO 6873, Type 4; e.g. Fujirock II, GC). Impressions can be copperor silver-plated.

ADA recommends following Council on Dental Materials, Instruments and Equipment. Disinfection of Impressions. JADA 1991;122(3):110; www.ada.org

Cartridge handling

The material is dispensed and mixed automatically by extruding through a mixing tip. Remove the transport seal before using a new cartridge. Then attach the mixing tip and, if required, an intraoral tip. Following use, the mixing tip should be used as a seal and should be left on the cartridge until the next application. The mixing tip should only be replaced with a new one when the cartridge is to be used again.

Recommended use

Monophase impression

Insert the impression tray loaded with Detail Mono into the patient's mouth. Maintain this position until the stipulated intraoral setting time of the material has elapsed. After removing the impression, thoroughly clean it and rinse under running water for approx. 15 s. Thoroughly dry the impression to obtain a bubblefree model.

For the Version in Maxi Dispenser of 380 ml

Insert Detail Mono NEXt cartridge in the mixing device. When doing so observe the user instructions issued by the device manufacturer.

1. Open the locking clip by moving it along the guide to its end position, and then remove the original cap or used mixing tip. (Figure 1-2)
2. Now position the new mixing tip and press both openings on the rear side of the mixing tip firmly on to the respective outlet openings in the cartridge. (Figure 3)
3. Once positioned, lock the mixing tip in place with the aid of the locking clip, by moving it along its guide until it surrounds the mixing tip and reaches the end position (Figure 4)

Note 1: Always fit a new mixing tip before each mixing procedure!

The mixing unit is unable to extrude impression material through set material in the mixing tip.

Note 2: Ensure that the mixing tip is fitted securely to the cartridge! Otherwise impression material will be pressed over the mixing tip and may cause contamination.

Technical data

• Dispenser Version (Fast Set)

Working time (incl. mixing time) ≤ 1:30 min
Recommended time in mouth ≥ 2:30 min
Strain in compression 2,2 – 3,5%
Compression set ≤ 0,5%
Shore-A hardness (after 24 h) ≤ 60

• Maxi Dispenser Version (Regular)

Working time (incl. mixing time) ≤ 2min
Recommended time in mouth ≥ 3:30 min
Mixing Ratio (base/catalyst) 5:1
Strain in compression 2,0 – 2,60%
Compression set ≤ 0,5%
Shore-A hardness (after 24 h) ≤ 65

Times refer to normal room temperature (23°C) and a normal relative humidity (50%). Higher temperatures shorten, lower temperatures increase those times.

Composition

Vinyl polysiloxanes (VPS), silicon dioxide, platinum catalyst, tensides and in small amounts stabilizers and additives.

Classification

DIN EN ISO 4823:2000, Type 2, Medium bodied consistency

Remarks

Store impressions at room temperature and normal humidity. Store Detail Mono NEXt at room temperature (15-25°C). Do not use after the expiry date.

Packaging

Detail Mono NEXt in Dispenser

2 cartridges of 50 ml each, 12 mixing tips, 12 intraoral tips

Detail Mono NEXt in Maxi Dispenser

1 maxi cartridge of 380 ml, 10 mixing tips.

For dental use only. Keep out of reach of children.

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Classification

DIN EN ISO 4823:2000, Type 2, consistance Medium bodied.

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