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INSTRUCTIONS FOR USE

General

AXIA® FILL N.F. is a light-curing radiopaque dental filling-material for all classes of cavities. It is based on Bis-GMA resins and inorganic fine glass filler particles in particle sizes of approximately 0.7 µm. The high filler content and the optimized resin matrix of AXIA® FILL N.F. ensure high compressive strength, very good color stability, excellent abrasion resistance, good polishability and very small polymerization shrinkage.

Composition

AXIA® FILL N.F. contains Barium glass (silan.) 74,2%, pyrog. silica (silan.) 2,8% in a matrix of Bis-GMA, TEDMA and UDDMA. Total inorganic filler content: 77 Wt. %.

Coloration

AXIA® FILL N.F. colors are based on the internationally recognized Vita® color system. Vita® is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co KG.

Contraindications

The application of composites is contraindicated in the following cases:

- when it is not possible to obtain a dry working area or to follow the foreseen application technique,
- for direct application on dentine near the pulp and when the pulp chamber has been exposed,
- if a possible allergy against one of the components of AXIA® FILL N.F. is known.

Side effects

Systemic side effects are not known. In some cases allergic reactions of type IV have been described.

Precautions for use

A suitable pulp/dentin protection is necessary to prevent pulp reactions. The work area should be kept dry. Contact with skin and eyes should be avoided. If accidental contact occurs, rinse immediately with plenty of water.

Interactions with other substances

Do not use AXIA® FILL N.F. together with materials containing eugenol or clove oil as they could inhibit the polymerization and cause discolorations. Water and oily air must be avoided, because they inhibit the polymerization at the contact area.

Cationic mouth washes as well as plaque-regulators and chlorohexidine may lead to discolorations.

Dosage and application

Isolation: The use of rubber-dam is recommended.

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www.dentalica.com – E-mail: dentalica@dentalica.com

Cleaning

Clean the tooth surface with a rubber cup and fluoride- and oil-free pumice paste. Remove residues completely with water spray and dry with oil-free air.

Shade selection

After cleaning, choose the color of the patient natural teeth with the Vita® shade guide. If possible, confront the color sample with the tooth at day light and while the tooth is still wet.

Cavity preparation

Prepare the cavity in the usual manner with minimal tooth reduction and retaining where possible the occlusal areas. Preparation of undercuts is not necessary unless great compressive or tensile forces act on the filling. To achieve lasting margins it is recommended to bevel the enamel edges, especially in the occlusal area.

Clean well with water spray and dry with oil-free air. The cavity should be dry and free of residues, however, do not desiccate dentin.

For class II cavities, interproximal wedges and matrix bands are indicated for a strong interproximal contact.

Pulp protection

To avoid possible pulp reactions in deep cavities cover the dentin nearest to the pulp with a thin layer of AXIA® PULP or other calcium hydroxide lining material.

Application of a liner

If a liner is used see instructions of the material.

Enamel etching

Apply AXIA® ETCH (etching gel) onto the enamel edges. Please follow the instructions of the product, then rinse and dry with oil and water-free air. It is essential that the etched enamel is not contaminated by saliva or any other substance until the procedure is carried on (e.g. until the adhesive is put).

Etching precaution

Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water. Refer to a physician if necessary!

Application of an adhesive

It is recommended to apply an adhesive (e.g. ONE-Q-BOND C.G.T. one-component adhesive) to the dry etched area. In this case see the separate instructions for use of the adhesive.

Application of AXIA® FILL N.F.

Apply the chosen paste directly into the cavity using a nonmetallic instrument. Light-cure in layers of 1.5 to 2 mm with an appropriate light unit to minimize polymerization shrinkage. Light-cure the first layer from buccal or lingual with a light-wedge. If a metal matrix is used the layer should only be 1 mm thick, because it can only be irradiated from the occlusal surface. Repeat this technique until the cavity is filled. Cure each layer for 30 to 40 seconds. When using matrix bands, do not loosen them during application or curing. If it is impossible to cure the composite under a matrix, the oxygen inhibited surface layer will be slightly softer. This layer must be trimmed off in the conventional manner.

Finishing

Contour with fine diamonds, stones or burrs. Polish with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips. With fine disks and silicone-polishers a high gloss surface may be obtained.

Classification

AXIA® FILL N.F. meets the requirements of international standards ISO 4049. Light-curing depth (40 seconds of lightening): no less than 3.5 mm.

Storage

The shelf-life of the paste is 3 years when stored at ≤ 25 °C (78 °F). Close immediately the syringe after use to avoid curing of the material by ambient light. Do not use after expiry date.

Packaging

AXIA® FILL N.F. Kit

6 syringes of 3 g each in the colors:

A1, A2, A3, A3.5, A4, B2, B3, C2, C3, IT, D.A2, D.A3, D.B3, D.C3.

AXIA® FILL N.F. Refills

1 syringe of 3 g in one of the following colors:

A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, C4, D3, IT, D.A2, D.A3, D.B3, D.C3, Bleach L, Beach XL.

For dental use only. Keep out of reach of children.

PRECISION products are carefully developed and extensively tested. We therefore guarantee their impeccable quality when used in accordance with their instructions for use. Application, use and processing of the product are beyond our control and are subject to the responsibility of the user.

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