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PRECISION

AXIA® FLUID

Light-curing flowable composite
for Class III, IV and V restorations

INSTRUCTIONS FOR USE

General

AXIA® FLUID is a low-viscosity light-curing flowable composite, containing fluoride, for Class III, IV and V restorations or as a liner in Class I and II cavities.
Thanks to its fluidity, AXIA® FLUID perfectly bonds to the cavity walls and reduces polymerization contraction stress. This leads to a smooth surface with optimal integrity.

Composition

AXIA® FLUID is made up of a Bis-GMA resin and other dimethacrylate resins, with a Barium glass inorganic charge and fluoride. It also contains catalysts, stabilizers and pigments.

Coloration

AXIA® FLUID colors are based on the internationally recognized Vita® color system. Vita® is a registered trademark of Vita Zahnfabrik H. Rauter GmbH & Co KG.

Contraindications

- The application of composites is contraindicated in the following cases:
- when it is not possible to obtain a dry field or follow the foreseen application technique,
 - for direct application on dentine near the pulp and when the pulp chamber is exposed,
 - if a possible allergy to the components of AXIA® FLUID is known.

Side effects

None known to date. In singular cases especially predisposed persons might develop hypersensitivity against single components of the material.

Precautions for use

To avoid possible pulp reactions cover the dentin nearest to the pulp with a thin layer of a suitable material. The working field must be dry (use of a rubber dam is recommended). Avoid the contact of the material with skin and eyes. In case of contact, rinse immediately with copious water.

Interactions with other substances

Do not use AXIA® FLUID with products containing eugenol or clove oil as they could inhibit the composite polymerization and cause changes in color. Avoid air spray containing oil or water as they could avoid the

composite polymerization in the contact points. Cationic collutories, clorexidine and plaque locators can cause changes in color.

APPLICATION METHOD

Isolation: Rubber dam is the recommended method of isolation.

Cleaning of elements

Before preparation, clean the tooth surface with a rubber cup, and fluoride- and fat-free pumice paste. Remove residues completely with water spray and dry with oil-free air.

Color choice

After cleaning, choose the color of the patient natural teeth with the Vita® shade guide. If possible, confront the color sample with the tooth at day light and while the tooth is still wet.

Cavity preparation

Prepare the cavity in the usual manner. Preparation of undercuts is not necessary in small cavities. Clean thoroughly the cavity with water spray and dry with oil-free air. The cavity should be dry and free of residues, however, do not desiccate the dentin.

Pulp protection

To avoid possible pulp reactions in deep cavities cover the dentin nearest to the pulp with a thin layer of AXIA® PULP or other calcium hydroxide lining material (follow attentively the producer instructions).

Enamel and dentin etching

AXIA® FLUID is compatible with all the common etching techniques. The application of AXIA® ETCH, etching gel with 37% phosphoric acid, is recommended starting from enamels surfaces (30 sec. at least) and continuing with the dentin ones (15 sec.). Then rinse with water for 30 seconds and dry with oil and water free air, without desiccating the dentin. The etched enamel should have a chalky white appearance. If this is not the case, etching has to be repeated until the desired effect is reached (read attentively the directions for use for the etching gel too).

Etching precaution

Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water. Réfer to a physician if necessary.

Application of adhesive

To enhance the adhesion of AXIA® FLUID the application of an adhesive (g.e. ONE-Q-BOND C.G.T. one-component adhesive) on the treated surface is recommended. In this case, follow attentively the indications for use.

Application of AXIA® FLUID

Apply AXIA® FLUID directly with the point of the syringe depositing a small drop in the base of the preparation and along the axial and gingival margins. During reconstruction use the syringe like a brush to spread AXIA® FLUID on the marginal surfaces. Light-cure the first quantity

applied for 30 seconds. For a perfect polymerization apply AXIA® FLUID in layers not thicker than 1-2 mm. Thinner layers should be light-cured for no more than 15-20 seconds. Shape the reconstruction margins using the same application point of the syringe for the maximum control. Finally light-cure for further 30 seconds.

Finishing and polishing

With AXIA® FLUID all the polishing standard techniques can be followed.

Classification

AXIA® FLUID meets the requirements of international standards ISO-4049.

Light-curing depth (40 seconds of lightening): 1-2 mm.

Storage

Shelf life is 3 years if stocked at a temperature not higher than 25°C. Do not use after expiry date.

Packaging

AXIA® FLUID

1 syringe of 3 g, available in the following colors: A1, A2, A3, A3.5, B3, C2, C3; 15 application tips.

For dental use only. Keep out of reach of children.

PRECISION products are carefully developed and extensively tested. We therefore guarantee their impeccable quality when used in accordance with their instructions for use. Application, use and processing of the product are beyond our control and are subject to the responsibility of the user.

For more information please contact:
DENTALICA S.p.A. – 22 via Rimini – I-20142 Milan – Italy
tel. #39-02895981 – fax #39-0289504249
email: dentalica@dentalica.com – http://www.dentalica.com

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DENTALICA S.p.A. – 22 via Rimini - I - 20142 Milan – Italie