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PRECISION

IMPLACEM

Non-eugenol,
acrylic-urethane polymer based cement
for implant luting

INSTRUCTIONS FOR USE

CE
0482

General

IMPLACEM is a non-eugenol cement designed in order to cement, permanently too, crown and bridges on implant abutments. IMPLACEM is unique because it is different from the other cements used for fixing implants. Most current non-eugenol temporary cements are compounds that contain zinc oxide and oil. These materials, while familiar, still have many limitations such as high film thickness, high solubility and difficulty in cleaning up the excesses and their progressive polymerization that reduces the possibility to remove the bridge after some time. IMPLACEM does not contain essential oils or zinc oxide, so it will not harm temporary crowns or interfere with permanent cementation. IMPLACEM is compatible with both acrylic and composite (bis-acryl) temporary materials and is indicated on both natural and metal implant abutments. In addition to being a versatile cement, IMPLACEM is ideal for permanent cementation of implant-retained crowns due to its polymerization stability and reduced solubility at the margins. IMPLACEM ensures easy cementation and high-retention and when it is properly used conserves, it maintains retrievability when desired.

Composition and characteristics

The essential components of IMPLACEM are: Bis-GMA, urethan dymethacrylate, multifunctional méthacrylates, polymerization activators. Solubility: less than 0.04 mg/cm³. Film thickness: less than 10 micron. Resistance to compression: 65 Mpa.

Indications

IMPLACEM is indicated for the permanent cementation of implant-retained crowns and bridges.

Method of application

a) Mixing of IMPLACEM

1. Test the automix syringe before mounting the mixing tip, extruding a little quantity of material.
2. Mount the mixing tip on the cartridge and before applying the material on the crown, test the syringe again.

b) Cementation of crown on implants

Before cementing a permanent crown in metal/metal-ceramic, on implant abutments, check fit, occlusion, proximal contacts and complete all finishing procedures.

1. Dry and isolate the prepared abutments. Absolute dryness isn't necessary, but obvious contamination should be avoided.
2. Apply on the external and internal surface an isolating medium (vaseline) which serves to reduce the strength of the adhesion to the porosity of the surface of the crown thus allowing a better remotion of the excess material.
In case of temporary resin crowns (or crowns that will later be removed) polish properly the surface of the crown or on the surface of the core and apply a thin layer of the isolating medium (vaseline) also to the internal surface of the crown.
3. If the perimetral edges are superior to 3.0 mm, it is helpful to place a retraction cord as per impression taking, in order to prevent that cement traces go in depth.
4. Insure that internal surfaces of the crown are clean and dry.
Apply a minimum quantity of the mixed cement, limited only to the edges of the crown in order to obtain a better control of the excess. Working time is 45-60 seconds.
5. Position immediately the crown, close the mouth until it reaches its centric occlusion. If it does not lift up, close it with a cotton roll waiting the initial cement set of 90-120 seconds. The initial set is characterized by a distinct rubbery phase in which excess cement can be easily removed using an explorer or periodontal probe.
6. Remove excess cement carefully, to take care to pass a probe through interproximal embrasures to eliminate excess cement in those areas, too.
7. Final cement set is between 4-5 minutes.
8. Properly control the absence of excess in the subgingival area.

Note

If cementing a temporary crown over a freshly placed composite resin core, we recommend very lightly lubricating the core material with a very thin film of a water soluble lubricant to prevent sticking. This step is not absolutely indicated, but certain brands of composite core material may have unpolymerized regions for several hours after curing which may cause cementation to IMPLACEM.

Storage

Store the product between 8°C and 20°C . Do not use after expiration date. Refrigerate for longer shelf life.

Packaging

2 automix syringues of 5ml each
15 tips for application

For dental use only. Keep out of reach of children.

PRECISION products are carefully developed and extensively tested. We therefore guarantee their impeccable quality when used in accordance with their instructions for use. Application, use and processing of the product are beyond our control and are subject to the responsibility of the user.

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DENTALICA S.p.A. – Via Rimini 22 – I-20142 Milano – Italia